

AVASURE GUARDIAN ACCESS DEVICE

AVA-B-0105

INSTALLATION INSTRUCTIONS

P/N: A-031701 REV. 1.0



COMPONENTS


Additional items may be required to comply with building code.



1 x Guardian Access Device (Part# AVA-B-0105)




1 x Wall Bracket (Part# A-026803)



2 x M3 x 10 mm T-10 6-Lobe
Pan Head Thread Forming Screws
(Part# A-012919)

TOOLS REQUIRED



TORX T-10 Bit or Screwdriver

SETUP

STEP 1

On wall, within a reasonable distance to hospital television, align provided metal bracket to appropriate height. Mount Wall Bracket a minimum of 63 in [1600 mm] from bottom of bracket to the floor.



STEP 2

Using proper screws, per substrate, attach bracket to wall. Use no less than a quantity of 3, 2 in [50.8 mm] minimum length, no less than #10 screws.

Do not use screws with flatheads. Screws must penetrate a wood or metal stud, .75 in [19 mm] thick plywood blocking, or anchors into poured concrete or concrete blocks.

STEP 3

Attach device to bracket by aligning plastic wire management housing, part of the pivot mechanism, to the slots located on the top and bottom flanges of wall bracket. The housing will drop into slotted top flange and snap into slotted bottom flange.



STEP 4

Secure device using supplied 10 mm long, M3, pan head, tamper-resistant, Torx drive screws. Pass screws through the 2 holes in bottom flange of wall bracket. Hand tighten screws into corresponding holes in bottom of housing.

STEP 5

Organize required cables similar to image shown. Attach zip-ties to the wire management housing as needed. Power device using either POE++ ethernet cable or power adapter AvaSure P/N A-034001. Connect device to TV using an HDMI cable (not supplied). Verify with the Observer Client software that the device is working properly.



DOCUMENTATION AND CONFIGURATION

To complete installation, contact AvaSure Support at support@avasure.com or **844.428.2797**, place a ticket with the admin, and provide the following for each room:

1. Hospital Name and Location (City/State)
2. Contact name/number
3. Room Number(s) from STEP 4
4. Serial Number(s) of the AvaSure AVA-B-0105
5. Camera Serial Number(s) from STEP 4

All room numbers must have the corresponding AvaSure serial number and camera serial number.

AvaSure Support may provide a list of cameras that may require additional troubleshooting as necessary. AvaSure Support will contact the installation team once the cameras have been configured and tested.

AvaSure